

Dipeptidyl peptidase 4

Catalog No: #AP82456



Package Size: #AP82456-1 50ug #AP82456-2 100ug #AP82456-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Dipeptidyl peptidase 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Dipeptidyl peptidase IV, DPP IV, T-cell activation antigen CD26, Thymocyte-activating molecule, THAM, CD26, Cd26
Accession No.	Uniprot ID:P28843
Uniprot	P28843
GeneID	13482;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Dpp4

Note: This product is for in vitro research use only