

Endothelin-1

Catalog No: #AP82457

Package Size: #AP82457-1 50ug #AP82457-2 100ug #AP82457-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Endothelin-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	ET-1,Preproendothelin-1,PPET1
Accession No.	Uniprot ID:P13206
Uniprot	P13206
GeneID	403424;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:EDN1

Note: This product is for in vitro research use only