

## Sphingolipid delta(4)-desaturase DES1

Catalog No: #AP82459



Package Size: #AP82459-1 50ug #AP82459-2 100ug #AP82459-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Sphingolipid delta(4)-desaturase DES1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Cell migration-inducing gene 15 protein, Degenerative spermatocyte homolog 1, Membrane lipid desaturase, DES1, MLD, MIG15
Accession No.	Uniprot ID:O15121
Uniprot	O15121
GeneID	8560;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: DEGS1

Note: This product is for in vitro research use only