

Deoxycytidine kinase

Catalog No: #AP82465



Package Size: #AP82465-1 50ug #AP82465-2 100ug #AP82465-3 1mg

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Deoxycytidine kinase |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Rt |
| Other Names | dCK |
| Accession No. | Uniprot ID:P48769 |
| Uniprot | P48769 |
| GeneID | 79127; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Dck

Note: This product is for in vitro research use only