

N(G),N(G)-dimethylarginine dimethylaminohydrolase 1

Catalog No: #AP82466



Package Size: #AP82466-1 50ug #AP82466-2 100ug #AP82466-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	N(G),N(G)-dimethylarginine dimethylaminohydrolase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	DDAH-1,DDAHI,Dimethylargininase-1
Accession No.	Uniprot ID:Q9CWS0
Uniprot	Q9CWS0
GeneID	69219;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Ddah1

Note: This product is for in vitro research use only