

ATP-dependent RNA helicase DDX42

Catalog No: #AP82467



Package Size: #AP82467-1 50ug #AP82467-2 100ug #AP82467-3 1mg

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	ATP-dependent RNA helicase DDX42
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	DEAD box protein 42,RCJMB04_2e15
Accession No.	Uniprot ID:Q5F485
Uniprot	Q5F485
GeneID	419959;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:DDX42

Note: This product is for in vitro research use only