

## Developmental pluripotency-associated protein 2

Catalog No: #AP82469



Package Size: #AP82469-1 50ug #AP82469-2 100ug #AP82469-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Developmental pluripotency-associated protein 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Pluripotent embryonic stem cell-related gene 1 protein,PESCRG1
Accession No.	Uniprot ID:Q7Z7J5
Uniprot	Q7Z7J5
GeneID	151871;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:DPPA2

Note: This product is for in vitro research use only