

Desmoglein-1

Catalog No: #AP82470



Package Size: #AP82470-1 50ug #AP82470-2 100ug #AP82470-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Desmoglein-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Cadherin family member 4,Desmosomal glycoprotein 1,DG1,Pemphigus foliaceus antigen,CDHF4
Accession No.	Uniprot ID:Q02413
Uniprot	Q02413
GeneID	1828;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:DSG1

Note: This product is for in vitro research use only