

Elongation factor G, mitochondrial

Catalog No: #AP82472



Package Size: #AP82472-1 50ug #AP82472-2 100ug #AP82472-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Elongation factor G, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	EF-Gmt,Elongation factor G 1, mitochondrial,mEF-G 1,Elongation factor G1,hEFG1,EFG,EFG1,GFM
Accession No.	Uniprot ID:Q96RP9
Uniprot	Q96RP9
GeneID	85476;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:GFM1

Note: This product is for in vitro research use only