

Sphingolipid delta(4)-desaturase DES1

Catalog No: #AP82473



Package Size: #AP82473-1 50ug #AP82473-2 100ug #AP82473-3 1mg

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Sphingolipid delta(4)-desaturase DES1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	Degenerative spermatocyte homolog 1,DES1,RCJMB04_22e12
Accession No.	Uniprot ID:Q5F3C1
Uniprot	Q5F3C1
GeneID	421327;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:DEGS1

Note: This product is for in vitro research use only