

Dachshund homolog 1

Catalog No: #AP82475



Package Size: #AP82475-1 50ug #AP82475-2 100ug #AP82475-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Dachshund homolog 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Dach1,Dach
Accession No.	Uniprot ID:Q9QYB2
Uniprot	Q9QYB2
GeneID	13134;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Dach1

Note: This product is for in vitro research use only