

N(G),N(G)-dimethylarginine dimethylaminohydrolase 2

Catalog No: #AP82477



Package Size: #AP82477-1 50ug #AP82477-2 100ug #AP82477-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	N(G),N(G)-dimethylarginine dimethylaminohydrolase 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	DDAH-2,DDAHII,Dimethylargininase-2,Protein G6a,S-phase protein,DDAH,G6A,NG30
Accession No.	Uniprot ID:O95865
Uniprot	O95865
GeneID	23564;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: DDAH2

Note: This product is for in vitro research use only