

PPAR α (Phospho-Ser21) Antibody

Catalog No: #12432

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	PPAR α (Phospho-Ser21) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography using non-phosphopeptide.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	PPAR α (Phospho-Ser21) Antibody detects endogenous levels of PPAR α only when phosphorylated at Ser21
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human PPAR α (Phospho-Ser21)
Target Name	PPAR α
Modification	Phospho
Other Names	PPARA, PPAR Alpha, PPARalpha, PPAR-alpha, NR1C1, HPPAR, PPAR
Accession No.	Swiss-Prot#: Q07869NCBI Gene ID: 5465
Uniprot	Q07869
GeneID	5465;
Target Species	human
Calculated MW	52kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000

Note: This product is for in vitro research use only