

Beta-defensin 4

Catalog No: #AP82480



Package Size: #AP82480-1 50ug #AP82480-2 100ug #AP82480-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Beta-defensin 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	BNBD-4,BNDB-4,BNBD4
Accession No.	Uniprot ID:P46162
Uniprot	P46162
GeneID	286836;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:DEFB4

Note: This product is for in vitro research use only