

Granulocyte colony-stimulating factor

Catalog No: #AP82481



Package Size: #AP82481-1 50ug #AP82481-2 100ug #AP82481-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Granulocyte colony-stimulating factor |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Pig |
| Other Names | G-CSF |
| Accession No. | Uniprot ID:O02837 |
| Uniprot | O02837 |
| GeneID | 396839; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: CSF3

Note: This product is for in vitro research use only