

Gamma-enolase

Catalog No: #AP82484



Package Size: #AP82484-1 50ug #AP82484-2 100ug #AP82484-3 1mg

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Gamma-enolase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	2-phospho-D-glycerate hydro-lyase,Enolase 2,Neural enolase,Neuron-specific enolase,NSE,Eno-2
Accession No.	Uniprot ID:P17183
Uniprot	P17183
GeneID	13807;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Eno2

Note: This product is for in vitro research use only