

Endothelial protein C receptor

Catalog No: #AP82485



Package Size: #AP82485-1 50ug #AP82485-2 100ug #AP82485-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Endothelial protein C receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Activated protein C receptor,APC receptor,Endothelial cell protein C receptor,CD201
Accession No.	Uniprot ID:Q4V811
Uniprot	Q4V811
GeneID	362248;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Procr

Note: This product is for in vitro research use only