

Sorbitol dehydrogenase

Catalog No: #AP82487



Package Size: #AP82487-1 50ug #AP82487-2 100ug #AP82487-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Sorbitol dehydrogenase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	L-iditol 2-dehydrogenase,Sdh1
Accession No.	Uniprot ID:Q64442
Uniprot	Q64442
GeneID	20322;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Sord

Note: This product is for in vitro research use only