

DNA fragmentation factor subunit beta

Catalog No: #AP82492



Package Size: #AP82492-1 50ug #AP82492-2 100ug #AP82492-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	DNA fragmentation factor subunit beta
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Caspase-activated deoxyribonuclease,CAD,Caspase-activated nuclease,CPAN,DNA fragmentation factor 40 kDa subunit,DFF-40,CAD,DFF2,DFF40
Accession No.	Uniprot ID:O76075
Uniprot	O76075
GeneID	1677;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:DFFB

Note: This product is for in vitro research use only