

## Cytoplasmic FMR1-interacting protein 2

Catalog No: #AP82493



Package Size: #AP82493-1 50ug #AP82493-2 100ug #AP82493-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Cytoplasmic FMR1-interacting protein 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	p53-inducible protein 121,KIAA1168,PIR121
Accession No.	Uniprot ID:Q96F07
Uniprot	Q96F07
GeneID	26999;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: CYFIP2

Note: This product is for in vitro research use only