

Epsin-3

Catalog No: #AP82494



Package Size: #AP82494-1 50ug #AP82494-2 100ug #AP82494-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Epsin-3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	EPS-15-interacting protein 3
Accession No.	Uniprot ID:Q9H201
Uniprot	Q9H201
GeneID	55040;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: EPN3

Note: This product is for in vitro research use only