

Discoidin domain-containing receptor 2

Catalog No: #AP82499



Package Size: #AP82499-1 50ug #AP82499-2 100ug #AP82499-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Discoidin domain-containing receptor 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Discoidin domain receptor 2,CD167 antigen-like family member B,Discoidin domain-containing receptor tyrosine kinase 2,Neurotrophic tyrosine kinase, receptor-related 3,Receptor protein-tyrosine kinase TKT,Tyrosine-protein kinase TYRO10,CD167b,NTRKR3,TKT,TYRO10
Accession No.	Uniprot ID:Q16832
Uniprot	Q16832
GeneID	4921;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:DDR2

Note: This product is for in vitro research use only