

Platelet basic protein

Catalog No: #AP82501



Package Size: #AP82501-1 50ug #AP82501-2 100ug #AP82501-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Platelet basic protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	PBP,C-X-C motif chemokine 7,Small-inducible cytokine B7,CXCL7,SCYB7
Accession No.	Uniprot ID:P43030
Uniprot	P43030
GeneID	396870;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PPBP

Note: This product is for in vitro research use only