

DnaJ homolog subfamily B member 1

Catalog No: #AP82504

Package Size: #AP82504-1 50ug #AP82504-2 100ug #AP82504-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	DnaJ homolog subfamily B member 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Heat shock 40 kDa protein 1,HSP40,Hsp40,Hspf1
Accession No.	Uniprot ID:Q9QYJ3
Uniprot	Q9QYJ3
GeneID	81489;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Dnajb1

Note: This product is for in vitro research use only