

Endothelin-1 receptor

Catalog No: #AP82507



Package Size: #AP82507-1 50ug #AP82507-2 100ug #AP82507-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Endothelin-1 receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Endothelin A receptor,ET-A
Accession No.	Uniprot ID:P26684
Uniprot	P26684
GeneID	24326;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ednra

Note: This product is for in vitro research use only