

Delta-like protein 4

Catalog No: #AP82508



Package Size: #AP82508-1 50ug #AP82508-2 100ug #AP82508-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Delta-like protein 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Drosophila Delta homolog 4,Delta4
Accession No.	Uniprot ID:Q9JI71
Uniprot	Q9JI71
GeneID	54485;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:DII4

Note: This product is for in vitro research use only