

Collagen alpha-1(XXVIII) chain

Catalog No: #AP82513

Package Size: #AP82513-1 50ug #AP82513-2 100ug #AP82513-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Collagen alpha-1(XXVIII) chain
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Col28
Accession No.	Uniprot ID:Q2UY11
Uniprot	Q2UY11
GeneID	213945;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Col28a1

Note: This product is for in vitro research use only