

Beta-defensin 4

Catalog No: #AP82514



Package Size: #AP82514-1 50ug #AP82514-2 100ug #AP82514-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Beta-defensin 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	BD-4,Defensin, beta 4,Bdef4
Accession No.	Uniprot ID:P82019
Uniprot	P82019
GeneID	56519;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Defb4

Note: This product is for in vitro research use only