

# Cullin-5

Catalog No: #AP82519



Package Size: #AP82519-1 50ug #AP82519-2 100ug #AP82519-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Cullin-5
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	CUL-5,Vasopressin-activated calcium-mobilizing receptor 1,VACM-1,VACM1
Accession No.	Uniprot ID:Q93034
Uniprot	Q93034
GeneID	8065;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:CUL5

Note: This product is for in vitro research use only