

DNA fragmentation factor subunit beta

Catalog No: #AP82521



Package Size: #AP82521-1 50ug #AP82521-2 100ug #AP82521-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	DNA fragmentation factor subunit beta
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Caspase-activated deoxyribonuclease, CAD, DNA fragmentation factor 40 kDa subunit, DFF-40, Cad
Accession No.	Uniprot ID:Q99N34
Uniprot	Q99N34
GeneID	84359;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Dffb

Note: This product is for in vitro research use only