

## Epoxide hydrolase 4

Catalog No: #AP82526



Package Size: #AP82526-1 50ug #AP82526-2 100ug #AP82526-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Epoxide hydrolase 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Abhydrolase domain-containing protein 7, Epoxide hydrolase-related protein, Abhd7, Ephxrp
Accession No.	Uniprot ID:Q6IE26
Uniprot	Q6IE26
GeneID	384214;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: Ephx4

Note: This product is for in vitro research use only