

# Interferon-induced, double-stranded RNA-activated protein kinase

Catalog No: #AP82529

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP82529-1 50ug #AP82529-2 100ug #AP82529-3 1mg

## Description

Product Name	Interferon-induced, double-stranded RNA-activated protein kinase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Eukaryotic translation initiation factor 2-alpha kinase 2,eIF-2A protein kinase 2,Interferon-inducible RNA-dependent protein kinase,Protein kinase RNA-activated,PKR,Tyrosine-protein kinase EIF2AK2,Prkr
Accession No.	Uniprot ID:Q63184
Uniprot	Q63184
GeneID	54287;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Eif2ak2

Note: This product is for in vitro research use only