

Delta(3,5)-Delta(2,4)-dienoyl-CoA isomerase, mitochondrial

Catalog No: #AP82536

Package Size: #AP82536-1 50ug #AP82536-2 100ug #AP82536-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Delta(3,5)-Delta(2,4)-dienoyl-CoA isomerase, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Accession No.	Uniprot ID:Q62651
Uniprot	Q62651
GeneID	64526;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Ech1

Note: This product is for in vitro research use only