

## Beta-defensin 1

Catalog No: #AP82538



Package Size: #AP82538-1 50ug #AP82538-2 100ug #AP82538-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |  |
|--------------------|--|
| Product Name       | Beta-defensin 1  |
| Brief Description  | Recombinant Protein                                    |
| Host Species       | E.coli   |
| Purification       | Greater than 90% by SDS-PAGE                           |
| Species Reactivity | Ms   |
| Other Names        | BD-1,Defensin, beta 1                                  |
| Accession No.      | Uniprot ID:P56386                                      |
| Uniprot            | P56386   |
| GeneID             | 13214;   |
| Tag Info           | His  |
| Storage            | Store at -20°C. (Avoid repeated freezing and thawing.) |

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Defb1

Note: This product is for in vitro research use only