

Diacylglycerol kinase alpha

Catalog No: #AP82546

Package Size: #AP82546-1 50ug #AP82546-2 100ug #AP82546-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Diacylglycerol kinase alpha
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	DAG kinase alpha, 80 kDa diacylglycerol kinase, Diglyceride kinase alpha, DGK-alpha, DAGK, DAGK1
Accession No.	Uniprot ID:P20192
Uniprot	P20192
GeneID	397097;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: DGKA

Note: This product is for in vitro research use only