

Dihydropteridine reductase

Catalog No: #AP82547

Package Size: #AP82547-1 50ug #AP82547-2 100ug #AP82547-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Dihydropteridine reductase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	HDHPR, Quinoid dihydropteridine reductase, Dhpr
Accession No.	Uniprot ID:Q8BVI4
Uniprot	Q8BVI4
GeneID	110391;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Qdpr

Note: This product is for in vitro research use only