

Ectonucleotide pyrophosphatase/phosphodiesterase family member 1



Catalog No: #AP82557

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP82557-1 50ug #AP82557-2 100ug #AP82557-3 1mg

Description

Product Name	Ectonucleotide pyrophosphatase/phosphodiesterase family member 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	E-NPP 1,Lymphocyte antigen 41,Ly-41,Phosphodiesterase I/nucleotide pyrophosphatase 1,Plasma-cell membrane glycoprotein PC-1,Npps,Pc1,Pdnp1
Accession No.	Uniprot ID:P06802
Uniprot	P06802
GeneID	18605;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Enpp1

Note: This product is for in vitro research use only