

DNA (cytosine-5)-methyltransferase 1

Catalog No: #AP82562



Package Size: #AP82562-1 50ug #AP82562-2 100ug #AP82562-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	DNA (cytosine-5)-methyltransferase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	Dnmt1, DNA methyltransferase Ggal, DNA MTase Ggal, MCMT, AIM, DNMT
Accession No.	Uniprot ID:Q92072
Uniprot	Q92072
GeneID	396011;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: DNMT1

Note: This product is for in vitro research use only