

Enoyl-CoA hydratase, mitochondrial

Catalog No: #AP82564



Package Size: #AP82564-1 50ug #AP82564-2 100ug #AP82564-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Enoyl-CoA hydratase, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Enoyl-CoA hydratase 1, Short-chain enoyl-CoA hydratase, SCEH
Accession No.	Uniprot ID:Q8BH95
Uniprot	Q8BH95
GeneID	93747;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Echs1

Note: This product is for in vitro research use only