

Diacylglycerol kinase alpha

Catalog No: #AP82567



Package Size: #AP82567-1 50ug #AP82567-2 100ug #AP82567-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Diacylglycerol kinase alpha
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	DAG kinase alpha,80 kDa diacylglycerol kinase,Diglyceride kinase alpha,DGK-alpha,Dagk,Dagk1
Accession No.	Uniprot ID:P51556
Uniprot	P51556
GeneID	140866;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Dgka

Note: This product is for in vitro research use only