

Fatty acid-binding protein, liver

Catalog No: #AP82572



Package Size: #AP82572-1 50ug #AP82572-2 100ug #AP82572-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Fatty acid-binding protein, liver
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	Fatty acid-binding protein 1,Liver basic FABP,LB-FABP,Liver bile acid-binding protein,L-BABP,Liver-type fatty acid-binding protein,L-FABP
Accession No.	Uniprot ID:P80226
Uniprot	P80226
GeneID	395345;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:FABP1

Note: This product is for in vitro research use only