

Dickkopf-related protein 3

Catalog No: #AP82574



Package Size: #AP82574-1 50ug #AP82574-2 100ug #AP82574-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Dickkopf-related protein 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Dickkopf-3,REIC,UNQ258/PRO295
Accession No.	Uniprot ID:Q9UBP4
Uniprot	Q9UBP4
GeneID	27122;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:DKK3

Note: This product is for in vitro research use only