

## C-X-C chemokine receptor type 4

Catalog No: #AP82576



Package Size: #AP82576-1 50ug #AP82576-2 100ug #AP82576-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	C-X-C chemokine receptor type 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	CXC-R4,Fusin,Leukocyte-derived seven transmembrane domain receptor,LESTR,Pre-B-cell-derived chemokine receptor,PB-CKR,Stromal cell-derived factor 1 receptor,SDF-1 receptor,CD184,Cmkar4,Lestr,Sdf1r
Accession No.	Uniprot ID:P70658
Uniprot	P70658
GeneID	12767;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Cxcr4

Note: This product is for in vitro research use only