

Testosterone 17-beta-dehydrogenase 3

Catalog No: #AP82577



Package Size: #AP82577-1 50ug #AP82577-2 100ug #AP82577-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Testosterone 17-beta-dehydrogenase 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	17-beta-hydroxysteroid dehydrogenase type 3,17-beta-HSD 3,Short chain dehydrogenase/reductase family 12C member 2,Testicular 17-beta-hydroxysteroid dehydrogenase,EDH17B3,SDR12C2
Accession No.	Uniprot ID:P37058
Uniprot	P37058
GeneID	3293;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:HSD17B3

Note: This product is for in vitro research use only