

Glutamate decarboxylase 2

Catalog No: #AP82582



Package Size: #AP82582-1 50ug #AP82582-2 100ug #AP82582-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Glutamate decarboxylase 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	65 kDa glutamic acid decarboxylase, GAD-65, Glutamate decarboxylase 65 kDa isoform, GAD65
Accession No.	Uniprot ID:Q05329
Uniprot	Q05329
GeneID	2572;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: GAD2

Note: This product is for in vitro research use only