

## Choline transporter-like protein 4

Catalog No: #AP82584



Package Size: #AP82584-1 50ug #AP82584-2 100ug #AP82584-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Choline transporter-like protein 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Solute carrier family 44 member 4,Ctl4,Ng22
Accession No.	Uniprot ID:Q91VA1
Uniprot	Q91VA1
GeneID	70129;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Slc44a4

Note: This product is for in vitro research use only