

Cystatin-A

Catalog No: #AP82589



Package Size: #AP82589-1 50ug #AP82589-2 100ug #AP82589-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Cystatin-A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Cystatin-A1, Stefin-A, Csta1
Accession No.	Uniprot ID:P56567
Uniprot	P56567
GeneID	209294;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Csta

Note: This product is for in vitro research use only