

Endoglin

Catalog No: #AP82591



Package Size: #AP82591-1 50ug #AP82591-2 100ug #AP82591-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Endoglin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	CD105
Accession No.	Uniprot ID:P37176
Uniprot	P37176
GeneID	397096;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ENG

Note: This product is for in vitro research use only