

Transcriptional repressor CTCF

Catalog No: #AP82592



Package Size: #AP82592-1 50ug #AP82592-2 100ug #AP82592-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Transcriptional repressor CTCF
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	11-zinc finger protein,CCCTC-binding factor,CTCF paralog
Accession No.	Uniprot ID:Q08705
Uniprot	Q08705
GeneID	396274;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CTCF

Note: This product is for in vitro research use only