

Interferon-induced, double-stranded RNA-activated protein kinase

Catalog No: #AP82609

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP82609-1 50ug #AP82609-2 100ug #AP82609-3 1mg

Description

Product Name	Interferon-induced, double-stranded RNA-activated protein kinase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Eukaryotic translation initiation factor 2-alpha kinase 2,eIF-2A protein kinase 2,Interferon-inducible RNA-dependent protein kinase,P1/eIF-2A protein kinase,Protein kinase RNA-activated,PKR,Serine/threonine-protein kinase TIK,Tyrosine-protein kinase EIF2AK2,p68 kinase,Pkr,Prkr,Tik
Accession No.	Uniprot ID:Q03963
Uniprot	Q03963
GeneID	19106;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Eif2ak2

Note: This product is for in vitro research use only